

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\TA_VB17\ChrisG\TRT_7T\1p5_TR3s_TE17ms_p3_6/8_axial_676b

TA: 15:15 PAT: 3 Voxel size: 1.5x1.5x1.5 mm Rel. SNR: 1.00 USER: ep2d_bold_sine_676b

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slice group 1	
Slices	70
Dist. factor	0 %
Position	L2.3 A38.5 H6.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.5 mm
TR	3000 ms
TE	17 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	A24

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	300
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
------------------	-------------

Series

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

V24	Off
A24	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Off

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L2.3 A38.5 H6.0
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	60 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	16
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Active
Meas[10]	Active
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

Sequence

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Introduction	Off
Bandwidth	1116 Hz/Px
Echo spacing	1 ms
<hr/>	
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
iPAT reference mode	Multi Shot
iPAT ref. averages	1
Readout type	Trapezoidal
Phase correction	Local
FatSat flip angle	110 deg
RF pulse type	Sinc
Excitation length	3072 us
BW-Time product	5.22
FFT factor	1.00